

Hybrid Kernelized Expectation Maximization implementation in STIR

Daniel Deidda^{1,2}, Ashley Gillman³

¹Division of Biomedical Imaging, School of Medicine, University of Leeds, UK.

²Department of Statistics, School of Mathematics, University of Leeds, UK.

³Commonwealth Scientific and Industrial Research Organisation, Brisbane, Australia.

Kernelised OSEM

$$\lambda_j = \sum_{f=1}^N \alpha_f k_{fj}^{(n)} \quad (1)$$

Iteratively reconstruct the coefficient α

$$\alpha_f^{(n+1)} = \frac{\alpha_f^{(n)}}{\sum_j k_{fj}^{(n)} \sum_i p_{ij}} \sum_j k_{fj}^{(n)} \sum_i p_{ij} \frac{1}{\sum_l p_{il} \sum_f k_{fl}^{(n)} \alpha_f^{(n)} + a_i} \quad (2)$$

The proposed kernel

$$k_{fj}^{(n)} = k_m(\mathbf{v}_f, \mathbf{v}_j) \cdot k_p(\mathbf{z}_f^{(n)}, \mathbf{z}_j^{(n)}) \quad (3)$$

$$k_m(\mathbf{v}_f, \mathbf{v}_j) = \exp\left(-\frac{\|\mathbf{v}_f - \mathbf{v}_j\|^2}{2\sigma_m^2}\right) \exp\left(-\frac{\|\mathbf{x}_f - \mathbf{x}_j\|^2}{2\sigma_{dm}^2}\right) \quad (4)$$

$$k_p(\mathbf{z}_f^{(n)}, \mathbf{z}_j^{(n)}) = \exp\left(-\frac{\|\mathbf{z}_f^{(n)} - \mathbf{z}_j^{(n)}\|^2}{2\sigma_p^2}\right) \exp\left(-\frac{\|\mathbf{x}_f - \mathbf{x}_j\|^2}{2\sigma_{dp}^2}\right) \quad (5)$$

New derived class and executable

- **KOSMAPOSL** derived from **OSMAPOSL**

Parameter file

KOSMAPOSLParameters:=

objective function type:= PoissonLogLikelihoodWithLinearModelForMeanAndList-ModeDataWithProjMatrixByBin

input:=

hybrid:= ;0 for KEM or 1 for HKEM

sigma_m:= ;default 1 controls anatomical edge preservation

sigma_p:= ;default 1 controls PET edge preservation

sigma_dm:= ;default: 1 (usual range 1-5) Spatial Gaussian scaling parameter for the anatomical prior (mm)

sigma_dp:= ;default: 1 (usual range 1-5) Spatial Gaussian scaling parameter for the PET prior (mm)

neighbours_num:= default:3 controls the size of the neighbourhood

num_non_zero_feat_elements:=default:1 controls the size of the feature vector patch

_image_filename:=if hybrid=1 and anatomical image is uniform the algorithm behave similarly to Snyder et al 1985

only_2D:=0

kernelised output filename prefix:=

Compatibility and functionality

- ▶ MR image has to be re-sliced to match the PET voxel size or vice versa (z voxel size depends on the scanner)
- ▶ The method would work for every scanner in STIR or customized (validated on Siemens mMR and mCT and GE signa).
- ▶ It is possible to reconstruct with both sinogram or list-mode, and use all the different objective functions
- ▶ If one wants to use any prior available within STIR this is also possible