

GPU Accelerated Time-of-Flight PET Reconstructions via STIR Integration of Parallelproj

N. Jurjew¹, D. Atkinson², K. Thielemans¹

1) Institute of Nuclear Medicine, UCL; 2) Centre for Medical Imaging, UCL

This research was supported by EPSRC Centre for Doctoral Training i4Health (EP/S021930/1) and Siemens Healthineers - Magnetic Resonance. Software used for this project is maintained with funding by CCP SyneRBI (EPSRC grant EP/T026693/1).

nicole.Jurjew.20@ucl.ac.uk

Introduction



Time of Flight (TOF) functionalities integrated into STIR^[1] master since v. 6.0

Traditional STIR Ray-Tracing projector does not support GPU utilization

Parallelproj^[2]: (independent) CUDA + C + Python library allowing use of CUDA GPUs for fast tomographic projection

Since STIR v. 5.0: wrapper for parallelproj projector into STIR for non-TOF reconstructions

This work

Extends the wrapper's functionality to support TOF data processing

[1] K. Thielemans et al., 'STIR: software for tomographic image reconstruction release 2', Physics in Medicine and Biology, Feb. 2012.

[2] G. Schramm and K. Thielemans, 'PARALLELPROJ—an open-source framework for fast calculation of projections in tomography', Front. Nucl. Med., Jan. 2024.

Materials & Methods (I)

PET phantom scans

Siemens Biograph Vision 600 [3, 4]

NEMA IEC Body phantom

Data generation – Siemens research software e7tools

Unlist LM data into sinograms

Estimate TOF scatter & randoms

Data generation – Prepare for use with STIR

Create STIR – readable headers (correct data order,...)

Use of Script to deal with this (STIR workshop at IEEE MIC 2023)

Conversion script to process e7tools output will be made available



[3] Casey, M.E., Osborne, D.R. (2020). Siemens Biograph Vision 600. In: Zhang, J., Knopp, M. (eds) Advances in PET. Springer, Cham.
[4] van Sluis J et al. Performance Characteristics of the Digital Biograph Vision PET/CT System. J Nucl Med. 2019 Jul; 60(7):1031-1036.

Materials & Methods (II)



Reconstructions

Using STIR via SIRF

MLEM, 10 iterations

For STIR Ray-Tracing: set number of LORs per bin to 5

Timing tools

Individual steps: STIR Timings Utility (stir_timings)

Full reconstructions: Using line-profiler ^[5]

A note on file sizes

prompts: 1.4 GB (signed integer)

randoms, scatter: 2.8 GB (single float)

Due to large file sizes: Caching switched off

[5] <https://pypi.org/project/line-profiler/>

System Specifications



CPU

- **Model:** Intel Xeon W2295
- **Cores/Threads:** 18 Cores / 36 Threads
- **Base/Boost Clock:** 3.0 GHz / 4.6 GHz

GPU

- **Model:** NVIDIA RTX A4000
- **CUDA Cores:** 6144
- **VRAM:** 16 GB GDDR6
- **Memory Bandwidth:** 448 GB/s
- **Base/Boost Clock:** 735 MHz / 1560 MHz

RAM

- **Size:** 64 GB DDR4 (2 x 32 GB modules)
- **Speed:** 2666 MT/s
- **Form Factor:** DIMM
- **ECC:** Supported

SW Versions

- **STIR:** 6.3
- **Parallelproj:** 1.9
- **OS:** Ubuntu 22.04.5 LTS (Jammy Jellyfish)
- **gcc:** 11.4.0
- **CUDA:** 12.1

A note on overhead

STIR: procedure for forward projection - copying data from CPU to GPU

- If image stored non-contiguously in memory, copy to contiguous memory block (rarely occurs!)
- Copy image from host to device (GPU)
- Call parallelproj (for all data, i.e. NOT subset-wise)
- Copy projection data from device to host
- Parallelproj and STIR projection data: different data order
 - Transpose TOF and spatial dimensions (on host) to internal ProjDataInMemory
- STIR is setup to forward project directly to a file
 - When data is accessed, it is copied from file into memory

Consequences:



- Double memory needed
 - Forward_project uses an internal buffer, but transposing uses its own buffer!
- Overhead in memory copies
- Overhead when only needing a subset of the data
 - can be avoided in SIRF (or STIR Python) by constructing one AcquisitionData per subset

Results – Timings



	RT [min]	PP [min]
forward projection	4.92	0.55
back projection	14.97	0.44
Full acquisition model set-up	16.67	1.50
Full reconstructions	133.33	5.00

Results – (STIR) Reconstructed Images


reconstructed image





reconstructed image coronal



Results – PP vs RT Forw. Projections




Simulated line-sources were
forward projected



Results – PP vs RT Forw. Projections

axial view




coronal view




Simulated line-sources were
forward projected


Relative Difference



Results – PP vs RT Forw. Projections




Simulated line-sources were
forward projected



Results – PP vs RT Forw. Projections


axial view



coronal view



Simulated line-sources were
forward projected



Acknowledgements




Special thanks go to Georg Schramm for his help with this work!

We would like to thank Judd Jones, Paul Schleyer and Mohammadreza

Teimoorisichani for their support!

Aside from valuable scientific input, they have been very supportive in making the e7tools conversion script publicly available!

This research was supported by EPSRC Centre for Doctoral Training i4Health (EP/S021930/1) and Siemens Healthineers - Magnetic Resonance. Software used for this project is maintained with funding by CCP SyneRBI (EPSRC grant EP/T026693/1).



References



- [1] K. Thielemans et al., 'STIR: software for tomographic image reconstruction release 2', Physics in Medicine and Biology, vol. 57, no. 4, pp. 867–883, Feb. 2012, doi: <https://doi.org/10.1088/0031-9155/57/4/867> .
- [2] G. Schramm and K. Thielemans, 'PARALLELPROJ—an open-source framework for fast calculation of projections in tomography', Front. Nucl. Med., vol. 3, Jan. 2024, doi: <https://doi.org/10.3389/fnume.2023.1324562>.
- [3] M. E. Casey and D. R. Osborne, 'Siemens Biograph Vision 600', in Advances in PET: The Latest in Instrumentation, Technology, and Clinical Practice, J. Zhang and M. V. Knopp, Eds., Cham: Springer International Publishing, 2020, pp. 71–91. doi: https://doi.org/10.1007/978-3-030-43040-5_6 .
- [4] J. van Sluis et al., 'Performance Characteristics of the Digital Biograph Vision PET/CT System', J Nucl Med, vol. 60, no. 7, pp. 1031–1036, Jul. 2019, doi: 10.2967/jnumed.118.215418.