

Respiratory motion correction of PET/MR data using SIRF

Paul S. Clegg, **Christoph Kolbitsch** and **Kris Thielemans**

School of Physics and Astronomy, University of Edinburgh

Physikalisch-Technische Bundesanstalt (PTB), Braunschweig and Berlin



University College London Institute of Nuclear Medicine



The problem

Einspänner *et al.* *EJNMMI Physics* (2022) 9:15

<https://zenodo.org/records/7862046>





Volumes
Positions
Activities

Are all wrong if you don't correct motion!

Data + Surrogate

- Here data from cardiovascular patients Siemens Biograph mMR
- 50 min PET scan with Ga-68 FAPI injected 35 min before scan
- Free breathing radial MRI scan lasting 4.1 min




- We require a “surrogate” for the position of the patient
- Here we take surrogates from the raw image data

PET and motion



- The PET data is quite noisy
- To extract motion:
 - Unpack listmode data into 0.5 second sinograms
 - Convert sinograms into low-resolution projections through the patient
 - Crop to include primarily the diaphragm
 - Use **principal components analysis** in time

Diaphragm
visible in raw
data




PET and motion

- Using 0.5 sec frames – surrogate is too noisy to use
- Instead use 1.5 sec frames three times
- Start at $t = 0, 0.5, 1.0$ sec – then create waited sum of **loadings**




PET and motion


Bin 1




Bin 3




Scatter and randoms estimated
from a portion of all data – not
for the bin



RTA




MRI and motion




Coil 3 has highest variance of the 18 coils

Value from centre of k-space is a good surrogate for the breathing motion – **PCA isn't required!**


MRI and motion



Displacement y-direction, slice: 19




Displacement y-direction, slice: 19





Displacement fields in the SI direction relative to bin 3.

Could use these transformations to warp PET bins

Displacement y-direction, slice: 19




MRI and motion








This is register, transform, add (RTA)

Motion corrected image reconstruction (MCIR)
gives improvement around heart

Noticeably sharper




Consistency





Arm edges good

Consistency

Surrogate comparison 2




Surrogate comparison 3



- Independent respiratory surrogates
- Agree perfectly without any adjustment

Breath-hold and the future



- Respiratory pattern far from regular
- Breath-hold protocols particularly challenging
- There is no StarVIBE MRI to accompany these segments

paul.clegg@ed.ac.uk

