


List-mode tracking and motion correction for PET imaging using low-activity fiducial markers

Marc Chamberland, MSc (PhD in Spring 2015!)

Robert deKemp, PhD


Tong Xu, PhD



STIR User's Meeting
at IEEE MIC



Carleton
UNIVERSITY



UNIVERSITY OF OTTAWA
HEART INSTITUTE
INSTITUT DE CARDIOLOGIE
DE L'UNIVERSITÉ D'OTTAWA

My claim to fame, from the STIR Wiki FAQ:


How do I add my own scanner to STIR?

You might not need to. If you specify the scanner geometry in your Interfile header, STIR will handle it ok.

For instance, you could use `create_projdata_template`, pick a scanner that might be somewhat similar to yours, and then edit the generated Interfile header. The scanner part of the header takes the same information as `Scanner::set_params()` (take care of changes between mm and cm). Obviously, it contains more information such as the actual number of views, ring differences etc that is supposed to be in your data. (Check the STIR Glossary as well for some info). Once you have this template, you should be good to go.

Alternatively, you will have to modify the `Scanner` class. Marc Chamberland gave a [good explanation of this on the stir-users list](#).

Note that STIR (at least 2.2 and earlier) ignores `view_offset` and block information.




For instance, you could use `create_projdata_template`, pick a scanner that might be somewhat similar to yours, and then edit the generated Interfile header. The scanner part of the header takes the same information as `Scanner::set_params()` (take care of changes between mm and cm). Obviously, it contains more information such as the actual number of views, ring differences etc that is supposed to be in your data. (Check the STIR Glossary as well for some info). Once you have this template, you should be good to go.

Alternatively, you will have to modify the `Scanner` class. Marc Chamberland gave a [good explanation of this on the stir-users list](#).

Note that STIR (at least 2.2 and earlier) ignores `view_offset` and block information.


Outline



Motion tracking and correction (Poster M10-5)




Use of STIR



Results

Motion tracking and correction

Positron emission tomography (PET) imaging is routinely used to assess myocardial perfusion...



Walls of the heart insert of a phantom, filled with rubidium-82 (^{82}Rb) tracer. (Coronal slice) Above: no motion.


...but respiratory, cardiac, and patient body motion create artefacts on the images.



Walls of the heart insert of a phantom, filled with rubidium-82 (^{82}Rb) tracer. (Coronal slice) Above: no motion; *below: motion blurring*.



A fiducial point-like positron-emitting marker can be localized by finding the position in space which minimizes the root mean square distance to its coincidence lines.




We developed an iterative expectation-maximization (EM) algorithm to track the 3D location of fiducial positron-emitting markers.

Some background rejection is needed to deal with the tracer activity in the patient's body. (Not shown in figure.)

EM Clustering Algorithm


COIL: COIncidence Lines




References

- T. Xu et al., Med. Phys. 33(7), 2598–2609 (2006).
M. Chamberland et al., Med. Phys. 38(2), 810–819 (2011).

We conducted a phantom study with a thorax phantom.




We transformed the phantom into a dynamic one!



**Phantom resting
on wooden board
with wheels**


**Motor with
rotating plate**

The marker was tracked using the raw list-mode data acquired during the PET scan.



The motion trace of the marker can be used to apply motion correction to the raw list-mode data.

- Reference marker position
- Marker position






Use of STIR

To check if the motion correction worked, I had two choices:

Use the GE Discovery 690's console to reconstruct the motion-corrected list-mode...



...or use STIR.



To reconstruct with STIR, I had to rely on projection data. Adding GE list-mode support was outside my C++ comfort zone.

make corresponding derived classes of `CListModeData`, `CListRecord` etc.
/ `InputFileFormat` class.



Notes for developers

If you want to add a new type of list mode data, you have to make corresponding derived classes of `CListModeData`, `CListRecord` etc. You also have to modify make sure that `read_from_file<CListModeData>` recognises your data. This normally involves creating a new `InputFileFormat` class.

Member Function Documentation

`virtual std::string stir::CListModeData::get_name () const`

pure virtual

Returns the name of the list mode data.


I added (non-validated) support for the GE Discovery 690 scanner in STIR by modifying Scanner.h...

Scanner.h

Scanner.h > class Scanner

```
92 class Scanner
93 {
94 public:
95
96     /** static members */
97     static Scanner * ask_parameters();
98
99     /** get the scanner pointer from the name
100     static Scanner * get_scanner_from_name(const string& name);
101     /** get the list of all names for the particular scanner
102     static string list_all_names();
103
104     // E931 HAS to be first, Unknown_scanner HAS to be last
105     // also, the list HAS to be consecutive (so DO NOT assign numbers here)
106     // finally, test_Scanner assumes E966 is the last in the list of CTI scanners
107     // supported by ecat_model from the LLN matrix library
108     // 08-3-2004, zlong, add user defined scanner
109     /** enum for all predefined scanners
110     /* \a Userdefined_scanner can be used to set arbitrary scanner geometry info.
111     \a Unknown_scanner will be used when parsing (e.g. from an Interfile header)
112     to flag up an error and do some guess work in trying to recognise the scanner from
113     any given parameters.
114     */
115     enum Type {E931, E951, E953, E921, E925, E961, E962, E966, E1080, Siemens_mMR, RPT, H1DAC,
116                Advance, DiscoveryLS, DiscoveryST, DiscoverySTE, DiscoveryRX, Discovery600, Discovery690,
117                HZLR, RATPET, PANDA, HYPERimage, nanoPET, HRRT, Allegro, GeminiTF, User_defined_scanner,
118                Unknown_scanner};
119
```


...and Scanner.cxx using the published characteristics of the scanner (and some detective work).




The image shows a screenshot of a C++ code editor with the file Scanner.cxx open. The code is part of a switch statement with two cases: Discovery690 and HZLR. Two red boxes highlight specific parts of the code. The first box, around lines 332-334, contains the expression $826.70F/2.F - 8.4F$. An arrow points from the text 'Effective ring diameter and average depth of interaction.' to this box. The second box, around lines 337-343, contains a list of integers: 4, 2, 6, 9, 1, 1, 1. An arrow points from the text 'Number of blocks per module, crystals per module, etc.' to this box.

```
324
325 case Discovery690:
326
327     set_params(Discovery690, string_list("GE Discovery 690", "Discovery 690", "GE D690", "D690"),
328               24,
329               381,
330               331, // TODO
331               2 * 288,
332               826.70F/2.F - 8.4F,
333               8.4F,
334               6.54F,
335               2.25446F,
336               static_cast<float>(-4.9730*_PI/180), //sign? TODO value
337               4,
338               2,
339               6,
340               9,
341               1,
342               1,
343               1);
344 break;
345
346 case HZLR:
347
348     set_params(HZLR, string_list("Positron HZL/R"),
349               32, 256, 2 * 192,
350               780.0F, 7.0F, 5.1875F, 2.F, 0.0F,
351               0, 0, 0, 0, 0, 0, 1);
352 // Default 7.0mm average interaction depth.
353 // crystals per singles unit etc unknown
354 break;
355
```

Effective ring diameter and average depth of interaction.

Number of blocks per module, crystals per module, etc.

The list_detector_and_bin_info utility was extremely useful in wrapping my head around STIR and projection data conventions.

A terminal window titled 'Marc — bash — 90x11' with standard macOS window controls (red, yellow, green buttons). The terminal shows the command 'list_detector_and_bin_info D690 250 500 12 17' and its output, which includes sampling distance, bin coordinates, LOR cylindrical coordinates, and detection position in Cartesian and index formats.

```

Last login: Wed Nov 12 11:15:14 on ttys000
[~] Marc$ list_detector_and_bin_info D690 250 500 12 17
sampling distance of adjacent LORs = 2.25446
bin: (segment 5, axial pos 12, view = 231, tangential_pos_num = 38)
bin coordinates: (tantheta: 0.0404199, m: 19.62, phi: 2.51982, s: 85.058)
LOR cylindrical: (z1: 3.26999, z2: 35.97, phi: 2.51982, beta: 0.207258 (= s: 85.058))
Detection position Cartesian: {3.26999, 378.344, 166.475}{35.97, -279.255, -304.754}
Detection position index (c:250,r:12,l:0)-(c:500,r:17,l:0)
[~] Marc$
```

This allowed me to check if the detectors were where I thought they were in Cartesian coordinates.


Fun discovery: everybody uses different conventions! (E.g. left-handed coordinate system instead of right-handed.)

In order to reconstruct the images, I wrote a utility that converts GE list-mode to projection data.

```
205 //assume it's a coincidence event
206 listmode.read((char*)&coincidenceEvent, sizeof(coincEvent));
207
208 //check if it's really a coincidence event
209 if (coincidenceEvent.coinc == 1) {
210     int X1 = coincidenceEvent.hiXtalRadialID;
211     int Z1 = coincidenceEvent.hiXtalAxialID;
212     int X2 = coincidenceEvent.loXtalRadialID;
213     int Z2 = coincidenceEvent.loXtalAxialID;
214
215     //in STIR, the first ring is at the back of the scanner;
216     //for GE, it is at the front
217     Z1 = (num_rings - 1) - Z1;
218     Z2 = (num_rings - 1) - Z2;
219
220     Bin bin;
221     proj_data_info_sptr->get_bin_for_det_pair(bin, X1, Z1, X2, Z2)
222
223     //for each coordinate, make sure it is positive
224     //i.e. segments and tangential bins can be negative in STIR
225     int segment = bin.segment_num() + (int)(num_segments/2);
226     int view = bin.view_num();
227     int z = bin.axial_pos_num();
228     int r = bin.tangential_pos_num() + (int)(num_tangential_pos/2);
229
230     //fill proj data array
231     ++projectionData4D[segment][view][z][r];
232 } //if (coinc event)
```

Member of
ProjDataInfoCylindricalNoArcCorr class

In order to reconstruct the images, I wrote a utility that converts GE list-mode to projection data.



I also did some work on a utility to convert crystal and geometric efficiency files from the D690 to projection data.

```
create_normalisation_files_from_GE_D690_calibration_files.cxx
main(int argc, char *argv[])

177 //now loop over projection data array and fill it with factors
178 for ( int segment = min_segment_num; segment <= max_segment_num; ++segment )
179 {
180     //each segment has a different number of axial coordinates
181     int num_axial_pos = proj_data_info_sptr->get_num_axial_pos(segment);
182     for ( int view = min_view_num; view <= max_view_num; ++view )
183     {
184         //a new array is needed every time because of the changing
185         //number of axial coordinates
186         array2D projectionData(boost::extents[num_axial_pos][num_tangential_pos]);
187         for ( int z = 0 ; z < num_axial_pos ; ++z)
188         {
189             for ( int r = min_tangential_pos_num ; r <= max_tangential_pos_num ; ++r)
190             {
191                 Bin bin = Bin(segment, view, z, r);
192                 int detector1, detector2;
193                 int ring1, ring2;
194
195                 proj_data_info_sptr->get_det_pair_for_bin(detector1, ring1, detector2, ring2, bin);
196
197                 //in STIR, the first ring is at the back of the scanner; for GE, it is at the front
198                 ring1 = (num_rings - 1) - ring1;
199                 ring2 = (num_rings - 1) - ring2;
200
201                 float detector1_efficiency = crystalEfficiency[ring1][detector1];
202                 float detector2_efficiency = crystalEfficiency[ring2][detector2];
203
204                 //in STIR, the tangential bin index can be negative, so I need to
205                 //shift it and make it positive to store it in the right element
206                 int r_index = r + (int)(num_tangential_pos/2);
207                 projectionData[z][r_index] = detector1_efficiency * detector2_efficiency;
```

**Member of
ProjDataInfoCylindricalNoArcCorr class**

Side note: I found neat stuff that is being worked on for future STIR releases.

file:///Users/Marc/STIR/STIR-doc/documentation/doxy/html/classstir_1_1LORAs2Points

Public Member Functions | List of all members


stir::LORAs2Points< coordT > Class Template Reference

Items related to Line Of Responses (preliminary)

A class for LORs. [More...](#)

```
#include "stir/LORCoordinates.h"
```

Inheritance diagram for stir::LORAs2Points< coordT >:



My tracking and correction code implements its own LOR class and methods...


Public Member Functions

```
const CartesianCoordinate3D
    < coordT > & p1 () const
CartesianCoordinate3D< coordT > & p1 ()
const CartesianCoordinate3D
    < coordT > & p2 () const
CartesianCoordinate3D< coordT > & p2 ()
LORAs2Points (const Ca
LORAs2Points (const LO
LORAs2Points (const LO
LORAs2Points (const LORInAxialAndNoArcCorrSinogramCoordinates< coordT > &)
virtual self_type * clone () const
virtual Succeeded change_representation (LORInCylinderCoordinates< coordT > &, const double radius) const
virtual Succeeded change_representation (LORInAxialAndNoArcCorrSinogramCoordinates< coordT > &, const double radius) const
virtual Succeeded change_representation (LORInAxialAndSinogramCoordinates< coordT > &, const double radius) const
virtual Succeeded get_intersections_with_cylinder (LORAs2Points< coordT > &, const double radius) const
```

I implemented a function to get the intersections between an N-sided polygon and a LOR...

Detailed Description


```
template<typename coordT>
class stir::LORAs2Points< coordT >
```




Results

**Qualitatively, the images are (relatively) motion free.
(Several corrections were not applied!)**


Static scan



Fast, irregular motion



**Fast, irregular motion;
motion corrected**



**Qualitatively, the images are (relatively) motion free.
(Several corrections were not applied!)**

No motion correction




0-15 s

105-120 s

210-225 s

285-300 s



With motion correction

NOTE: Those images were reconstructed from the GE console, not by STIR.

Take-away message:

- 1. Motion tracking using low-activity fiducial positron-emitting markers and correcting the raw list-mode data from a GE D690 scanner.**
- 2. Using STIR to convert the list-mode and calibration files to projection data, correct the data, and reconstruct for qualitative assessment of motion.**
- 3. It works! (And I have patient data, but not quite ready to show.)**

Future work:

- 1. Better corrections after applying motion correction.**
- 2. LOR-as-2-points-based reconstruction?**
- 3. Contribute something to STIR!**

Acknowledgements

1. **Dr. Rob de Kemp (University of Ottawa Heart Institute)**
2. **All the PET staff involved (University of Ottawa Heart Institute)**
3. **Dr. Charles Stearns for technical discussion of the D690 (GE Healthcare)**
4. **Dr. Tong Xu (Department of Physics, Carleton University)**
5. **IEEE NSS/MIC Conference Trainee Grant for travel support**
6. **Dr. Thielemans and Dr. Tsoumpas for having answered some of my questions on the STIR mailing list over the years.**

**Thank you
for your time!**

