

Connecting STIR and the SIMIND Monte Carlo Simulator

R. Gillen^{1,2}, S. Porter¹, M. Ljungberg³, K. Thielemans¹

¹ Institute of Nuclear Medicine, University College London Hospital, London, UK,

² Marine, Medical and Nuclear Department, National Physical Laboratory, Teddington, UK

³ Lund University, Medical Radiation Physics, Lund, Sweden

Acknowledgements: RG is completing a PhD with funding from the Industrial Strategy Challenge Fund via NPL. STIR is maintained with the support from the Computational Collaborative Project in Synergistic

Reconstruction for Biomedical Imaging, CCP SynerBI, UK EPSRC grant EP/T026693/1.


Thanks to Dr Kjell Erlandsson and Ashley Morahan for assistance in getting started with SIMIND.

Contact Rebecca.Gillen.18@ucl.ac.uk

Introduction - SIMIND

- Monte Carlo (MC) simulation software ¹
- Simulates typical SPECT systems only
 - Less configurable but faster to set up and run than GATE


```
Simulation started.: 2022:03:08 15:03:29
Simulation stopped.: 2022:03:08 15:05:17
Elapsed time.....: 0 h, 1 m and 48 s
DetectorHits.....: 1692645
DetectorHits/CPUsec: 15861
```



[1] M. Ljungberg and S.-E. Strand, "A Monte Carlo program for the simulation of scintillation camera characteristics," *Comput. Methods Programs Biomed.*, vol. 29, pp. 257–272, 1989.

Using SIMIND with STIR


- SIMIND can define emission and density data with pre-set shapes
- More flexibility possible using “voxel based” input datasets (e.g. from STIR)
- There are differences in conventions between STIR and SIMIND:
 - Coordinate systems, attenuation/density maps, accepted input data types, interfile headers...etc



SIMIND demo – input data and other required files

- Emission data: generate in STIR and save *.v file as *.smi
- Density data: generate in STIR and save *.v file as *.dmi
- Simulation parameters: *.smc file
- Other optional input files – see SIMIND manual


16 bit integer data



Command Prompt

C:\simind>

SIMIND demo – setting simulation parameters



Note: In this version of **simind**, the protective cover works only for photons that are impinging on the lower side of the detector.

Index-14: Integer flag indicating the type of phantom that is to be simulated. Available standard phantoms are

- 2 = Rectangular phantom,
- 3 = Vertical cylindrical phantom, and
- 4 = Horizontal cylindrical phantom.

Note: The cylindrical phantom can be elliptical if either Index-5 and -6 or Index-6 and -7 (depending on the direction of the phantom) are chosen not to be equal.

If the flag value is negative, the phantom is simulated using binary density maps. The file name storing the particular maps are defined in Option 14 of the main menu and should have the following extensions:

- 1 = An density distribution map (16-bit integers, 1000*density, extension: *.dmi)
- 2 = The Zubal Voxel-Man density distribution (8-bit, extension: *.dat)
- 3 = The Zubal brain density distribution (8-bit, extension: *.dat)
- 4 = The Zubal whole-body Voxel-Man density distribution (8-bit, extension: *.dat)
- 5 = Code-based Zubal phantom (8-bit, extension: *.dat)
- 6 = Code-based NCAT/XCAT phantom (32-bit float, extension: *.bin)
- 7 = Generic NCAT/XCAT phantom (32-bit float, extension: *.bin)
- 8 = Generic MCAT phantom (32-bit float, extension: *.bin)

The non-homogeneous phantom is further defined by Index-31 to -38, which are described below.

Note: Because the coordinate system for MCAT/NCAT axially is rotated by 180 degrees

c:\ Command Prompt - change simind Ug example.smc

```
15 - Source type.....: -1.0000
←[22;07H←[K      Index number.....:
←[2J←[H
┌────────────────────────────────────────────────────────────────────────────────┐
│ C H A N G E: Scintillation camera parameters                               │
└────────────────────────────────────────────────────────────────────────────────┘
16 - Shift source in x-direction.....cm:      0.0000
17 - Shift source in y-direction.....cm:      0.0000
18 - Shift source in z-direction.....cm:      0.0000
19 - Photon direction.....deg:      2.0000
20 - Upper window threshold.....keV: -15.0000
21 - Lower window threshold.....keV: -15.0000
22 - Energy resolution ...[140 keV]..... %:      8.8000
23 - Intrinsic resolution [140 keV].....cm:      0.3500
24 - Emitted photons per decay.....:      0.8430
25 - Source activity.....MBq:      1.0000
26 - Number of photon histories * 1E6.....:      10.0000
27 - keV/channel.....keV:      0.5000
28 - Pixel size in simulated image.....cm:      0.4000
29 - SPECT: No of projections.....:      120.0000
30 - SPECT: Rotation [0=-360,1=-180,2=360,3=180] :      0.0000
←[22;07H←[K      Index number.....:
```


rebecca@DESKTOP-HKF3GNS: /mnt/c/Users/Rebecca/Documents/Code/inm-WIPs/people/RebeccaGillen/SIRF_SIMIND/UG_examp... — □ ×

rebecca@DESKTOP-HKF3GNS: /mnt/c/Users/Rebecca/Documents/Code/inm-WIPs/people/RebeccaGillen/SIRF_SIMIND/UG_examples\$

I

SIMIND demo – running the simulation

C:\simind\simind.exe

```
rebecca@DESKTOP-HKF3GNS:/mnt/c/Users/Rebecca/Documents/Code/inm-WIPs/people/RebeccaGillen/SIRF_SIMIND/UG_examples$  
simind.exe simind_ug_example.smc output_filename/FD:density_data/FS:emission_data/NN:10/PX:0.4
```

SIMIND Monte Carlo Simulation Program V7.0

```
-----  
Phantom(S): h2o      Crystal...: nai      InputFile.: simind_ug_example  
Phantom(B): bone     BackScatt.: pmt      OutputFile: output_filename  
Collimator: pb_sb2   SourceRout: smap     SourceFile: emission_data  
Cover.....: al      ScoreRout.: scattwin  DensityMap: density_data  
-----
```


```
PhotonEnergy.....: 140      ge-legp  PhotonsPerProj....: 2139780  
EnergyResolution...: 8.8      SPECT   Activity.....: 1  
MaxScatterOrder....: 3       Xrays   DetectorLenght....: 19.35  
DetectorWidth.....: 27.65    BScatt  DetectorHeight....: 0.95  
UpperEnergyWindowTresh: 150.5 Random Distance to det...: 26  
LowerEnergyWindowTresh: 125.5 Photon ShiftSource X.....: 0  
PixelSize J.....: 0.4       Resolt  ShiftSource Y.....: 0  
PixelSize S.....: 0.4       Header  ShiftSource Z.....: 0  
HalfLength S.....: 19       SaveMap HalfLength P.....: 19  
Halfwidth S.....: 0.1       Halfwidth P.....: 11  
HalfHeight S.....: 0.1      HalfHeight P.....: 11  
SourceType.....: Integer2Map PhantomType.....: Integer2Map  
-----
```

```
-----  
Projection      Hits      Scatt/Tot      Progress  
118 of 120      545812      0.000      97 %  
-----
```

4 minutes later...

SIMIND output


- SIMIND output files:
 - .h00 [header file(s)]
 - Need to convert to STIR-readable form
 - .a00 [sinogram data file(s)]
 - .res [results file – summary of params]
- Other.... see SIMIND manual




Reconstructing SIMIND output with STIR

- Convert .h00 file to .hs file using script
- Run reconstruction as normal
- Remember to use STIR-specific attenuation map if AC applied

STIR simulation & STIR recon



SIMIND simulation & STIR recon



Examples and more resources:

- **See STIR SIMIND Tutorial Repository:**
 - https://github.com/samdporter/STIR_SIMIND
- **Test datasets** include example *.smc file for SIMIND, and *.par and .sh scripts to generate data, convert headers to STIR-friendly format, and do recons with STIR
- **Jupyter Notebook** for those who work with STIR via Python
- **Crib sheet** includes more detail on data types, commands, and recommendations for setting SIMIND parameters for configuring the simulation to work with STIR

Conclusions/Discussion

- It would be useful for a number of applications to combine STIR and SIMIND software systems
- We have established basic requirements to connect STIR & SIMIND
 - Can now use SIMIND perform MC simulations using STIR-generated voxelwise input data
 - Can also use STIR to reconstruct data simulated with SIMIND
- This work is ongoing and we welcome discussion or feedback

**Thank you for your
attention!**