

Out of the Lab and Into the Market

Preliminary Evaluation of STIR for Commercial SPECT Reconstruction Software

Spencer Manwell, Ph.D.
Imaging Physicist and Software Developer
Convergent Imaging Solutions

Daniel Deidda, Ph.D.
Senior Scientist, Nuclear Medicine Imaging
National Physical Laboratory, UK

Who we are

- We're a **Software as a Medical Device (SaMD)** manufacturer based out of Ottawa, Canada.
 - Operating since ~2007
 - Our primary product (**UniSyn Molecular Imaging**) is a radiology and nuclear medicine workstation used for reading and processing diagnostic images.
 - Commercial partnerships/integration with several PACS vendors and SaMD manufacturers.



Why STIR?

- Expanding our software portfolio to include SPECT image reconstruction.
 - Clinics with older scanners no longer receiving software support/upgrades
 - Clinics with scanners from multiple vendors looking for a harmonized solution.
- We don't have in-house expertise needed to start from scratch so we turned to open-source solutions.
- Considered:
 - **PyTomography**
 - **CASToR**
 - **STIR**



Why STIR?

- Our customers are largely non-academic image centers.
 - Interested in **standard functionality**, that is **reliable** and **high-performing**.
 - Require I/O in DICOM format.
- For SPECT recon, PyTomography, CASToR and STIR have similar capability.
 - PyTomography and STIR seem to have the most complete solution
 - recon algos,
 - data correction,
 - support for DICOM format
- STIR Advantage
 - Code maturity (stability)
 - Active user group (You!)



Challenges (Opportunities?) With STIR

- We requirement for projection data and reconstructed images in DICOM format.
 - Limited support for parsing DICOM projection data
 - No support for outputting reconstructed images in DICOM format.
- Support for attenuation corrected SPECT reconstruction.
 - STIR requires mu-map the same spacing/extents as the SPECT image.
 - Such mu-maps are not provided by most scanners / CT images rarely have the same extents/spacing as the SPECT image.
 - Mu-map generation is laborious.
- Support for modern 3D resolution recovery techniques.
 - SPECTUB utility supports 2D and 3D resolution recovery.
 - 2D approach introduces banding artifacts
 - 3D approach more accurate but can be very slow for some scanners/collimators.
 - Gaussian approach not well suited for high-energy emitters (I-131, Y-90)

Current Status



Conversion of DICOM Projection Data to Interfile

See ***SPECT_dicom_to_interfile*** utility (B. Thomas, D. Deidda & K. Theilemans).

Limited support for multi-head SPECT system (See Issues: #1102, #1184).

Various attributes were calculated corresponding only to one-detector head

- *num_of_projections*
- *rotation_radius*
- *extent_of_rotation*
- These issues were addressed by modifying the DICOM read to loop over all items (detector heads) in the *DetectorInformationSequence* and extracting the information for all heads.



Conversion of DICOM Projection Data to Interfile

Failed to parse DICOM data written using *Transfer Syntax* that included compression (Issue #1474).


- Issue was a non-starter for DICOM dataset produced using this syntax as the utility would fail.
- Modified the DICOM file read to detect compressed files, rewrite them using a non-compression transfer syntax and re-read.



Converting Interfile Images to DICOM

STIR has no existing utilities to produce a DICOM image following reconstruction.

- We're developing a new utility ***SPECT_interfile_to_DICOM.cxx*** to handle this task.
- Some information still missing (Issue #1436):
 - The create time of the reconstructed image
 - Series Description
 - Image ID



Example Recon Image 2D vs 3D Resolution Recovery


Tc-99m MDP Bone Scan

From a Siemens Symbia SPECT/CT

- OSEM 16i 16s, CTAC, TEW, RR
- Increased noise seen on the coronal view but not on the transverse for 2D vs 3D RR.

2D Resolution Recovery

3D Resolution Recovery



Looking ahead...

- Wrap up our work towards completing DICOM I/O support with STIR
 - Pull requests are coming – sorry Kris!
- Streamline the process of mu-map creation
 - Hopefully ease the STIR restrictions on matching image extents/spacing
- Investigate acceleration techniques for pre-computing projection matrix, especially in the presence of 3D resolution recovery.
- Build a UI component to wrap the reconstruction code.
- Validation in collaboration with our users.





Thank you!

Contact Info:

Spencer Manwell, Ph.D.

Imaging Physicist and Software Developer

spencer.manwell@convergentimaging.com